SoftSurgicalStrip™ Chirurgische Tupfer (DE) MD

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SoftSurgicalStrip™ Surgical Sponges (EN) MD

ATTENTION!
The product is to be used by qualified users only!
Use only with undamaged packaging!

Product Description

SoftSurgicalStrip™ is made of a soft, pliable, moisture-absorbing material that is used to absorb fluid, protect tissue, and reduce bleeding.
SoftSurgicalStrip™ is equipped with a safety thread to facilitate retrieval.
SoftSurgicalStrip™ is optionally available with X-ray contrast thread to allow localization of the swabs.

Purpose / Indication for use

SoftSurgicalStrip™ is used in surgical procedures in neurosurgery.
Its main functions are to remove fluid accumulation in surgical wounds, protect tissue from damage, and moisturize tissue during prolonged neurosurgical procedures. The swab is not intended to remain in the body after the wound has closed.

Contraindications

SoftSurgicalStrip™ is not intended to remain in the body after closure of the wound.
No other contraindications are known.

Target population

The intended patients are individuals undergoing neurosurgical intervention on the central nervous system. There are no further restrictions on the patient population regarding age, gender, ethnicity, etc.

Intended User

The user group consists of specialists trained in surgical procedures in neurosurgery.

Clinical Claims

Protection against tissue damage
Absorption of fluids (absorption capacity of $\geq 500\%$ of its own weight)
Maintaining tissue moisture during prolonged surgical procedures

Recommended technique / application

Removal: Once the package is opened, the swabs may become contaminated. Therefore, remove the swabs from the package only immediately before use and dispose of opened packages after the operation. Immediately after opening the package, check the completeness of the swabs according to the number indicated on the package.

Placement: After opening the surgical area, open the appropriate size SoftSurgicalStrip™ pack and moisten the swabs using sterile water or physiological solution. Swabs then placed according to standard techniques and requirements. Apply sterile water or physiological solution to swabs during the procedure to improve conformity to tissue contours and provide additional moisture to the tissue.

Count card: After removing the swabs from the surgical area, collect them by hooking the safety thread in the side notches on the count card. Re-check for completeness before closing the surgical area to ensure that all used swabs have been removed after surgical procedures.

Safety thread: The safety thread on the swabs is used only to locate the swabs and must not be used to pull them out.

X-ray contrast strip (optional): The X-ray contrast thread enables the swab to be localized in radiological imaging.

Warnings and safety instructions:

- The product is to be used by qualified personnel only.
- Product for ONE TIME USE - DO NOT reuse, repackage or re-sterilize the swabs. Used swabs must be considered potentially infectious medical waste after use. They must not be reused after application in order to prevent infections.
- DO NOT use swabs from damaged packaging.
- CHECK PACKAGING: Not sterile if damaged, wet, seal broken, or past expiration date.
- DO NOT leave the swab in the body.
- DO NOT apply swab dry to wound area.
- DO NOT allow swabs to dry out during application.
- DO NOT remove swabs from the package until immediately prior to use and discard opened packages after surgery to avoid contamination.
- Use the counting card or use alternative counting methods to ensure that all used swabs have been removed after surgical procedures.
- Do NOT use the safety thread on the swab to pull it out. Otherwise tissue damage may occur.
- The product can be used in combination with all standard neurosurgical instruments. Caution: Simultaneous use of the swabs with laser systems may result in ignition or sparking of the material, severing of the safety thread, and damage to the swab.

Disposal

Used swabs must be considered potentially infectious medical material after use. Disposal must be carried out appropriately in accordance with applicable hospital and national waste regulations.

Packaging / single use

SoftSurgicalStrip™ is supplied sterile unless this package has been damaged or opened.
Each package contains 10 pieces.
SoftSurgicalStrip™ sponges are to be used only once.

Shelf-life

The swabs must not be used after the expiry date (hourglass symbol) indicated on the packaging. Opened packages must no longer be used. Store in the original packaging.

STERILE R

The swabs have been radiation sterilized in the final packaging.

Storage and transport Instructions

SoftSurgicalStrip™ must be stored and transported in an environment that maintains packaging and purity.
- Protect from moisture
- Protect from light
- Store in a temperature range from 4°C to 50°C

SoftSurgicalStrip™ is manufactured by:

neuro-competence GmbH
Birkenhügel 5
D-02689 Sohland a.d. Spree, Germany
Tel +49 (0) 35936-33 99-0
Fax +49 (0) 35936-33 99-11
E-mail sohland@neuro-competence.de

Should a serious incident or undesirable side effects occur during the handling and use of the product that are not listed in these instructions for use, please report them to us with a detailed description of the incident or side effects.

Please use the contact details above for this purpose.

We would like to point out that all serious incidents related to the product must be reported not only to the manufacturer but also to the competent authority of the Member State in which the incident occurred (https://health.ec.europa.eu/medical-devices-sector/new-regulations/contacts_en).

Summary of Safety and Clinical Performance Report (SSCP):

In accordance with Article 32 of Regulation (EU) 2017/745 on Medical Devices, the manufacturer shall provide a Summary of Safety and Clinical Performance Report (SSCP) for this device. This report contains a summary of the clinical data, the performance evaluation, and information on the safety and benefit-risk balance of the device. The summary is publicly available via the EUDAMED database and can be viewed at the following link: <https://ec.europa.eu/tools/eudamed>

Please note that access to the respective product information is via the UDI-DI or the trade name of the device.

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